

Figure 4 Comparison of viscosity measurement for 1000 cP oil at 25 °C using three geometry setups: cone and plate, 30 mm; parallel plate, 30 mm; and bob and cup, 25 mm.




Figure 6 Viscosity as a function of torque for three Newtonian standard oils: 12,500; 1000; and 10 cP at 25 °C using bob and cup, 25 mm.




Figure 5 Three independent runs for 1000 cP oil at 25 °C using bob and cup, 25 mm.




Figure 7 Viscosity as a function of shear rate for 1000 cP oil at various temperatures: 10; 25; 50; and 80 °C using bob and cup, 25 mm.

Steady shear experiments were performed using the Windows™-based MICRA user interface software (Microsoft, Redmond, WA) (Figure 3) from high to low shear rate.

The resulting steady-state viscosity (Pa s) is plotted versus torque (mNm) in Figure 4 (torque is plotted on the abscissa rather than shear rate [s^{-1}] to normalize the data for the three experiments). As expected, all three measuring systems show similar data. The ability to use different sample geometries permits the user to expand the capacity and range of the instrument and/or accommodate a wider range of materials and sample sizes/amounts. The bob and cup system is well suited for low-

viscosity materials and/or low-shear-rate measurements since it provides more sample contact surface area and thus more torque signal. Three repeat runs are shown in Figure 5, where viscosity (Pa s) is plotted as a function of torque (mNm) for the bob and cup system, demonstrating the performance of the MERLIN viscometer over the working torque range of the instrument.

The bob and cup geometry was utilized to test several different viscosity standards covering the range 10–10,000 cP at 25 °C. All data are in acceptable ranges of viscosity, as shown in Figure 6. Again, the data are plotted versus torque to normalize the abscissa.

Viscosity (Pa s) at several set temperatures for the nominal 1000 cP oil is shown in Figure 7 as a function of shear rate (s^{-1}). The ability to control and/or systematically change the sample test temperature is paramount in making viable rheology and viscosity measure-

Table 2 Temperature response of the MERLIN viscometer with bob and cup, 25 mm

Instrument setpoint	Meter reading
10 °C	10.1 °C
25 °C	25.1 °C
50 °C	50.0 °C
80 °C	79.9 °C




Figure 8 Temperature profiles for heating and cooling between 10 and 80 °C at a rate of 1 °C/min. The sample is 1000 cP oil using bob and cup, 25 mm, at a shear rate of 10 s^{-1} .




Figure 9 Viscosity as a function of $1/T$ for heating between 10 and 80 °C at a rate of 1 °C/min. The sample is 1000 cP oil using bob and cup, 25 mm, at a shear rate of 10 s^{-1} .

ments. The viscometer's built-in Peltier temperature control system supports all measuring systems, including cone and plate, parallel plate, and bob and cup. The temperature accuracy of the instrument using the bob and cup measuring system is summarized in Table 2. It is clear that the temperature was well calibrated throughout the range investigated. A thermal enclosure is available to enhance the performance at high temperatures.

Another advantage of the MERLIN viscometer over other commercial viscometers is the ability to perform user-defined, controlled temperature ramps. This is vitally important for temperature-sensitive materials. The

activation energy E_a can be estimated from such an experiment. Figure 8 shows the temperature profiles for heating and cooling experiments between 10 and 80 °C, both at a rate of 1 °C/min. A bob and cup geometry was used for this test on 1000 cP oil at a constant shear rate of 10 sec^{-1} . It is clear that the temperature can be controlled well with the Peltier temperature control system. Figure 9 shows the viscosity as a function of $1/T$ for the heating experiment, which can be fit by the Arrhenius model, $\eta = C e^{k/T}$, where the constant k is related to the activation energy by E_a , with R the Boltzmann constant. From this experiment, the activation energy obtained from the heating




Figure 10 Shear stress as a function of time for a commercial ketchup at 0.1 s^{-1} using bob and cup, 25 mm, at 25 °C.




Figure 11 Viscosity as a function of shear rate for a commercial ketchup using bob and cup, 25 mm, at 25 °C.

run was 16.9 kJ/mol. Other controlled temperature ramps from 0.1 to over 2 °C/min are also available.

The MERLIN viscometer is also able to measure the dynamic yield stress or yield point of a sample, and the MICRA software supports this experimental mode of testing. Using the 25-mm bob and cup, a commercial ketchup sample was evaluated. Figure 10 shows the results plotted as shear stress (Pa) versus time (s). The yield stress experiment was performed by applying a 0.1 s^{-1} shear rate and measuring the resulting stress response. The maximum stress

value measured is the yield stress of the sample, approximately 14.7 Pa for this ketchup. This yield stress measuring capability is a significant improvement over other commercially available viscometers, which only provide viscosity at a single rotational speed. The MICRA software also allows for user-defined shear rate sweeps, where the shear rate is increased/decreased stepwise, permitting a complete flow curve profile to be determined. *Figure 11* shows the viscosity versus shear rate for the commercial ketchup sample in the shear rate range 1–1000 s⁻¹. The data reveal that the sample is highly shear

thinning (pseudoplastic), and clearly cannot be defined by a single viscosity value as obtained from other commercial viscometers.

Conclusion

This article demonstrates that the MERLIN viscometer has the capacity to easily generate precise data in terms of viscosity, flow curves, and yield stress. Its built-in Peltier temperature control system and ability to conduct temperature isothermal and ramp experiments distinguishes it from other viscometers. In addition, the fact that it is equipped standard with different

measuring systems enables the viscometer to accurately and routinely measure the viscosity of a wide range of samples under varied conditions.

Dr. Dao is Senior Staff Engineer, and Dr. Ye is Staff Engineer, REOLOGICA Instruments, 231 Crosswicks Rd., Bordentown, NJ 08505, U.S.A.; tel.: 609-298-2522; fax: 609-298-2795; e-mail: info@reologicainstruments.com. Mr. Hutchison is President, RheoSys, Plainsboro, NJ, U.S.A. Dr. Hedman is VP Software Development, REOLOGICA Instruments, Lund, Sweden.